

### Features & benefits

#### EMCCD Technology

Ultimate in Sensitivity from EMCCD gain – even single photon signals are amplified above the noise floor. Full QE of CCD chip is harnessed (no intensifier).

#### TE cooling to -95°C

Unparalleled elimination of EM-amplified darkcurrent noise.

#### RealGain™

Absolute EMCCD gain selectable directly from a linear and quantitative scale.

#### iCam

Unique innovation that empowers the EMCCD to operate with market-leading acquisition efficiency through live cell microscopy software.

#### > 65% QE from virtual phase sensor

Highly efficient photon collection. One window design.

#### Extended red response

Significantly higher sensitivity to red-emitting dyes such as CY5, mCherry, dsRed and Alexa680. Bose Einstein Condensation in NIR.

#### Fast speed readout

35MHz readout speed delivers 31 frames/sec at full megapixel resolution; 60 frames/sec when 2x2 binned.

#### UltraVac™\*1

Critical for sustained vacuum integrity and to maintain unequalled cooling and QE performance, year after year.

#### 8 x 8 µm pixel size (fully binnable)

Excellent balance of NyQuist resolution and photon collection.

#### Minimal Clock-Induced Charge

Unique pixel clocking parameters, yielding minimized spurious noise floor.

#### Enhanced Baseline Clamp

Essential for quantitative accuracy of dynamic measurements.

#### Negligible EM Gain ageing

No requirement for gain recalibration.

#### Built-in C-mount compatible shutter (optional)

Easy means to record control dark images – excellent for optimization of experimental set-up.

#### Cropped sensor mode

Specialised acquisition mode for continuous imaging with fast temporal resolution.

### “ High resolution Megapixel EMCCD ”

A proud member of Andor's iXon<sup>EM</sup>+ EMCCD range, the 885 benefits from the unique innovations and high-end performance specifications that have characterized this camera family as the industry's leading high-performance EMCCD.

The megapixel sensor format and 8 x 8µm pixel size of the 885 presents an attractive combination of field of view and resolution, offering excellent Nyquist over-sampling for cell microscopy.




When more light is available from the sample, the EMCCD gain can be completely switched off and the camera operated as a 'traditional' CCD camera. However if the camera is used in low-light conditions, the EMCCD gain can be applied to render it single photon sensitive, while maintaining a full resolution frame rate of 31 frames/sec. The absolute EM gain multiplication can be varied linearly from unity up to a thousand times directly via RealGain™, a true quantitative EM gain scale. Extended red QE response is ideally matched to popular red-emitting fluorophores and for imaging of Bose Einstein Condensates using NIR probe laser

#### Camera overview

Active Pixels	1004 x 1002
Pixel Size (W x H; µm)	8 x 8
Image Area (mm)	8 x 8
Active Area Pixel Well Depth (e <sup>-</sup> , typical)	30000
Gain Register pixel well depth (e <sup>-</sup> , typical)	80000
Max Readout Rate (MHz)	35
Frame Rate (frames per sec)	
Full resolution	31.4
2 x 2 binning	60.5
Read Noise (e <sup>-</sup> )	
@ 35 MHz	25
with EM Gain	<1

#### Quantum efficiency\*2



#### Peak Quantum Efficiency

CCD Type	Typical
VP @ 600nm	65%


## Technical specifications

### System characteristics

<b>Pixel Readout Rate (MHz)</b>	35, 27, 13
<b>Linearity (% maximum)<sup>*3</sup></b>	1
<b>Vertical Clock Speed (µs)</b>	0.5 to 1.9 (variable)
<b>Electron Multiplier Gain</b>	1 - 1000 times (software controlled) via RealGain™ control with temperature compensation
<b>Digitization</b>	14-bit @ 35, 27 & 13 MHz readout rate
<b>Dark Current (e/pix/sec)<sup>*4</sup></b>	
@ -70°C	0.028
@ -85°C	0.012
@ -95°C	0.005

### System Readout Noise (e<sup>-</sup>)<sup>\*5</sup>

	Typical	With Electron Multiplication
<b>35MHz through EMCCD amplifier</b>	25	<1
<b>27MHz through EMCCD amplifier</b>	22	<1
<b>13MHz through EMCCD amplifier</b>	12	<1



### Minimum sensor temperatures (typical)<sup>\*6</sup>

	DV option	DU option
<b>Air cooled (ambient air at 20°C)</b>	-70°C	-80°C
<b>Re-circulator (XW-RECR); (ambient air @ 20°C)</b>	-80°C	-90°C
<b>Water-cooled (@ 10 °C, 0.75 l / min)</b>	-85°C	-95°C


### Operating & storage conditions

<b>Operating Temperature</b>	0°C to 30°C ambient
<b>Relative Humidity</b>	< 70% (non-condensing)
<b>Storage Temperature</b>	-25°C to 55°C

### Power requirements<sup>\*7</sup>

- 1A @ +12V
- 0.3A @ -12V
- 3.0A @ +5V

### Full frame rate<sup>\*8</sup>



### Max frames per second<sup>\*9</sup>

Binning	Array size			
	1002 x 1002 (Full Frame)	512 x 512	256 x 256	128 x 128
<b>1 x 1</b>	31.4	60.3	113	213.7
<b>2 x 2</b>	60.5	113.6	201.6	367.7
<b>4 x 4</b>	112.6	204.1	331.0	568.2

### Computer requirements

#### Minimum

- Windows 2000 or XP operating system
- PCI 2.2 or PCI-X 1.0 compatible computer (PCI slot must have bus master capability)
- Available auxiliary internal power connector
- 25 MBytes free hard disc space

#### Recommended

- 3.2 GHz Pentium (or better) + 1 GB RAM
- SATA RAID 0 hard disc, e.g. Seagate Barracuda, Western Digital Caviar RE or Raptor, etc.


**NOTE:** In all cases the operating system should be on a separate hard drive and the hardware controller should be on a separate PCI bus

### Need more information? Please contact us at:

International Office	US Office
<b>Phone:</b> +44 28 9023 7126	<b>Phone:</b> 800.296.1579
<b>Fax:</b> +44 28 9031 0792	<b>Fax:</b> 860.290.9566
Japanese Office	Chinese Office
<b>Phone:</b> +81 3 3511 0659	<b>Phone:</b> +86-10-5129-4977
<b>Fax:</b> +81 3 3511 0662	<b>Fax:</b> +86-10-6445-5401

## Ordering information &amp; notes

To order the camera you require, please use the following ordering system:



E.g. a DV-885K-C00-#VP is an iXon<sup>EM</sup> + 885 camera with 3-stage peltier vacuum cooling and without internal shutter.

The iXon<sup>EM</sup> + 885 also requires one of the following software options:

<b>Andor Solis (i)</b>	A ready-to-run Windows 2000 or XP-based package with rich functionality for data acquisition and processing.
<b>Andor SDK</b>	A DLL driver and software development kit that let you create your own applications for the Andor Camera. Available for Windows 2000 or XP and Linux.
<b>Andor iQ</b>	A comprehensive multi-dimensional imaging software package. Offers tight synchronization of EMCCD with a comprehensive range of microscopy hardware, along with comprehensive rendering and analysis functionality. Modular architecture for best price/performance package on the market.
<b>Third party software compatibility</b>	Drivers are available so that the iXon <sup>EM</sup> + range can be operated through a large variety of third party imaging packages.

The following accessories are available for use with the iXon<sup>EM</sup> + 885:

<b>XW-RECR</b>	Re-circulator for enhanced cooling performance
<b>XW-CHIL-150</b>	Chiller/re-circulator for maximum cooling performance




Image of Human Umbilical Vein Endothelial Cells (HUVEC) taken with an iXon<sup>EM</sup> + 885.



Courtesy of Christophe Dehio and Matthias Truttmann, University of Basel.

Specifications are subject to change without notice

- ◆1 Assembled in a state-of-the-art Class 10,000 cleanroom facility, Andor's UltraVac™ vacuum process combines a permanent hermetic vacuum seal (no o-rings), with a stringent protocol to minimize outgassing, including use of proprietary materials. Outgassing is the release of trapped gases that would otherwise prove highly problematic for high-vacuum systems.
- ◆2 Quantum efficiency of the CCD sensor as measured by the CCD Manufacturer.
- ◆3 Linearity is measured from a plot of Counts vs. Signal up to the saturation point of the system. Linearity is expressed as a percentage deviation from a straight line fit.
- ◆4 The dark current measurement is averaged over the CCD area excluding any regions of blemishes.
- ◆5 System Readout noise is for the entire system. It is a combination of CCD readout noise and A/D noise. Measurement is for Single Pixel readout with the CCD at a temperature of -75°C and minimum exposure time under dark conditions. Under Electron Multiplying conditions, the effective system readout noise is reduced to sub 1e<sup>-</sup> levels. Noise values will change with pre-amplifier gain (PAG) selection. Values quoted are measured with highest available PAG setting.
- ◆6 The iXon<sup>EM</sup> + 885 can be ordered either with 3-stage or 4-stage peltier vacuum cooling, the 'DV' or 'DU' options respectively.
- ◆7 These power requirements are the maximum load that will be drawn from the computer for the camera head and controller card combined.
- ◆8 The graph shows the full frame rates possible when reading out the sensor at 35, 27 and 13 MHz pixel readout rates, and using 0.5µs per row vertical clock speed.
- ◆9 The max frames / second for iXon<sup>EM</sup> + imaging CCDs is the maximum speed at which the device can acquire images in a standard system. Shown are the frame rates at 35 MHz digitization rates for a range of binning or array size combinations. Measurements are shown for 0.5µs per row vertical clock speed. It also assumes internal trigger mode of operation.

## Dimensions

Weight: 3.1 kg [7 lb 1 oz]



## Notes:

1. The clearance from the C-mount face plate to the shutter is 6mm. Please ensure that when fitting a lens, to a system with a built in shutter, that it does not extend into the housing by more than 5mm.
2. There are mounting holes (1/4-20UNC) located on three sides of the camera. They are positioned centrally at a distance of 40mm from the front of the front face.

## Connections

