

Little Genius



- ▶ 25×0.2ml, 16×0.5ml and 24×0.2ml (8-strips) three blocks for your choice.
- ▶ 20×2 character LCD display software interface, PCR document editing, browsing and streamline the monitoring operation easy and clear.
- ▶ Document Link function and temperature, time function automatically modified to meet the complex set of procedures.
- ▶ Compact, light weight, economical and practical.



- TC-25/H
- TC-16/H
- TC-24/H

8-strip / 24-microplate

Little Genius		
TC-25/H	TC-16/H	TC-24/H
25 × 0.2ml	16 × 0.5ml	24 × 0.2ml 8-strip 24-microplate
4°C ~ 99°C		
≥ 3.0°C/sec		≥ 3.0°C/sec
≥ 2.0°C/sec		≥ 2.8°C/sec
≤ ± 0.4°C		
≤ ± 0.4°C	≤ ± 0.5°C	≤ ± 0.3°C
105°C		
—		
—		
20 × 2 LCD		
—		
99 files		
—		
9		
99		
AC220V 50Hz 200VA		
285 × 218 × 178 (L × W × H)		
3.2kg		
RS232		

Distributore:

Via di Grotte Portella 6/8
00044 Frascati RM
Tel. 06 72632108
Fax 06 72676527